



## OPTIPHARM Co., Ltd.

63, Osongsaengmyeong 6-ro, Osong-eup, Heungdeok-gu,  
Cheongju-si, Chungcheongbuk-do, 28158 Rep. of KOREA  
Tel: +82-43-249-7500, Fax: +82-43-249-7501  
www.optipharm.co.kr


---

# PAM(Porcine Alveolar Macrophage) Cells

Porcine alveolar macrophages (PAM) are useful for studying viral infections and immune response in pigs  
Optipharm's SPF miniature pigs were also tested negative for PRRSV, PCV2, CSFV real time PCR

Product : **Optipharm SPF Miniature Pig – PAM Cells**

Type/volume : **1x10<sup>7</sup> cells / 1.0 ml, 10 vials**



---

### Specifications

Description : Porcine Alveolar Macrophage Cells are derived from the lung of a euthanized 3~4 weeks old SPF Miniature pig.

Suitable for research areas related to virus pathogenesis and host cell interactions

Species : SPF Miniature Pig (Yucatan)

Tissue/Organ : Lung

Donor Age : 3~4 weeks old

Cell Morphology : Macrophage-like

Cell Type : 1x10<sup>7</sup> cells / 1.0 ml

Applications : For Research Use Only

Preservation : Storage Temperature: Liquid nitrogen vapor phase.

Handling : This material is sold for in-vitro use only in manufacturing and research and is not suitable for human use. It is the responsibility of the user to undertake sufficient verification and testing to determine the suitability of each product's application.